



Table 1: Dimensions

Mass Spectrometer System	Pole diam	A	B	C	X	R1	R2	R3	R4	MF
HALO 100 RC	6	138	211	161.4	35	36	67	59	58	DN 35-CF
HAL 200 RC	6	188	211	161.4	35	36	67	59	58	DN 35-CF
HALO 201 RC	6	205	211	177.6	35	36	67	59	58	DN 35-CF
HAL 201 RC	6	248	225	220.3	35	36	67	59	58	DN 35-CF
HALO 301 RC	6	205	225	177.6	35	70	34	204	139	DN 35-CF
HAL 301 RC	6	248	225	220.3	35	70	34	204	139	DN 35-CF

MSIU- Mass Spectrometer Interface Unit : 19 inch rack mountable unit.

Height: 88mm (2U) Width: 448mm Depth: 405mm + 76mm for connector clearance

Dimensions are in mm.

Lengths identified include clearance for gauge and filament wiring, and are minimum clearance dimensions.